## Claims:

1. A method for the non-invasive early detection of colon cancer and/or intestinal cancer precursor cells by means of mutational analysis of the genes for APC, K-ras,  $\beta$ -catenin and B-raf in a sample, characterized in that

the method comprises the following steps:

- collecting a stool and/or tissue sample,
- homogenizing the sample,
- obtaining DNA from the sample,
- performing an amplification reaction in the genes for APC, K-ras,  $\beta$ -catenin and B-raf, using the primers
  - s1 TTGCAGTTATGGTCAATACCC
  - as1 GTGCTCTCAGTATAAACAGGATAAG
  - s2 CCTCAAAAGGCTGCCACTTG
  - as2 CTGTGACACTGCTGGAACTTCGC
  - s3 AGCACCCTAGAACCAAATCCAGCAG
  - as3 TGGCATGGTTTGTCCAGGGC
  - s4 ACAAACCATGCCACCAAGCAGA
  - as4 GAGCACTCAGGCTGGATGAACAAG
  - s5 TTCCAGATGCTGATACTTTA
  - as5 CTGAATCATCTAATAGGTCC

for APC, the primers

- s CTGGTGGAGTATTTGATAGTG
- as TCTATTGTTGGATCATATTC

for K-ras, the primers

- s CTGATTTGATGGAGTTGGAC
- as CTTGAGTGAAGGACTGAGA

for  $\beta$ -catenin, and the primers

- s TGTATCACCATCTCCATATC
- as GCATTCTGATGACTTCTGGT

for B-raf,

wherein amplification products are formed, and

- performing a mutational analysis in the amplification products.
- 2. The method according to claim 1, characterized in that the detection of mutations in selected sections of the genes for APC, K-ras,  $\beta$ -catenin and B-raf is effected by means of a DNA chip, said DNA chip including probes for APC, K-ras,  $\beta$ -catenin and B-raf from those regions of the above-mentioned genes that are flanked by the primer sequences specified in claim 1.
- 3. The method according to claim 1 or 2, characterized in that the APC, K-ras,  $\beta$ -catenin and B-raf genes are accumulated from total DNA by hybridizing sequence-specific biotinylated oligonucleotides with the genes for APC, K-ras,  $\beta$ -catenin and B-raf using coupling of the biotin residue to streptavidin and subsequent separation via magnetic particles.
- 4. The method according to claims 1 to 3, characterized in that amplification products, especially PCR products, are separated in an agarose gel for control purposes prior to purification.
- 5. The method according to any of claims 1 to 4, characterized in that the mutational analysis of the PCR products is effected using electrophoretic techniques, preferably SSCP, alternatively by means of a chromatographic procedure, preferably an HPLC-based procedure.

- 6. The method according to the preceding claim, characterized in that detected mutagenic conformations of a single strand are isolated and optionally sequenced.
- 7. Primer sequences selected from the group comprising: the primers
  - s1 TTGCAGTTATGGTCAATACCC
  - as1 GTGCTCTCAGTATAAACAGGATAAG
  - s2 CCTCAAAAGGCTGCCACTTG
  - as2 CTGTGACACTGCTGGAACTTCGC
  - s3 AGCACCCTAGAACCAAATCCAGCAG
  - as3 TGGCATGGTTTGTCCAGGGC
  - s4 ACAAACCATGCCACCAAGCAGA
  - as4 GAGCACTCAGGCTGGATGAACAAG
  - s5 TTCCAGATGCTGATACTTTA
  - as5 CTGAATCATCTAATAGGTCC
  - or alternatively
  - s2 GAATCAGCTCCATCCAAGT
  - as2 TTTCTGCTATTTGCAGGGT
  - for APC, the primers
  - s CTGGTGGAGTATTTGATAGTG
  - as TCTATTGTTGGATCATATTCG
  - for K-ras, the primers
  - s CTGATTTGATGGAGTTGGAC
  - as CTTGAGTGAAGGACTGAGAA
  - for  $\beta$ -catenin, and the primers
  - s TGTATCACCATCTCCATATC
  - as GCATTCTGATGACTTCTGGT
  - for B-raf.
- 8. Use of the primer sequences according to claim 7 in mutational analysis, especially in the analysis of the APC, K-ras,  $\beta$ -catenin and B-raf genes.

9. A kit, comprising primers selected from the group comprising:

the primers

- s1 TTGCAGTTATGGTCAATACCC
- as1 GTGCTCTCAGTATAAACAGGATAAG
- s2 CCTCAAAAGGCTGCCACTTG
- as2 CTGTGACACTGCTGGAACTTCGC
- s3 AGCACCCTAGAACCAAATCCAGCAG
- as3 TGGCATGGTTTGTCCAGGGC
- s4 ACAAACCATGCCACCAAGCAGA
- as4 GAGCACTCAGGCTGGATGAACAAG
- s5 TTCCAGATGCTGATACTTTA
- as5 CTGAATCATCTAATAGGTCC
- or alternatively
- s2 GAATCAGCTCCATCCAAGT
- as2 TTTCTGCTATTTGCAGGGT
- for APC, the primers
- s · CTGGTGGAGTATTTGATAGTG
- as TCTATTGTTGGATCATATTCG

for K-ras, the primers

- s CTGATTTGATGGAGTTGGAC
- as CTTGAGTGAAGGACTGAGAA

for  $\beta$ -catenin, and the primers

- s TGTATCACCATCTCCATATC
- as GCATTCTGATGACTTCTGGT

for B-raf,

and optionally information relating to combining the contents of the kit.

10. Use of the kit according to claim 9 in the detection of colon cancer and/or colon cancer precursor cells.